

# QUICK. EFFICIENT. EFFECTIVE.

POLICE LINE DO NOT CROSS

**QuikClot®  
COMBAT  
GAUZE<sup>LE</sup>**

## IT'S QUIKCLOT® OR IT'S NOT

### *Experience Improves Performance*

Z-Medica® training kits and moulage trainer simulate real-world bleeding challenges, so users get hands-on experience. Get to know your QuikClot® Law Enforcement products, register at [Z-Medica.com/training](http://Z-Medica.com/training)



#273



#384



#529



#350



#261



#534



#539

### *Proprietary Technology*

QuikClot Combat Gauze® LE is a soft non-woven material impregnated with kaolin, and is available in a 3"x4 yds z-fold, packed in an easy-tear, vacuum sealed pouch, contains an X-ray detectable strip for easy identification, and has a five-year shelf life.


# THE FIVE LARGEST U.S. POLICE DEPARTMENTS AND ALL FIVE BRANCHES OF THE MILITARY TRUST QUIKCLLOT TO SAVE LIVES



## Proven safe and effective

Independent test results published in several respected journals demonstrate that QuikClot has consistent hemostatic achievement, reduction in blood loss and strong blood clot formation.<sup>2,3,4,5,6</sup>

## Consistent hemostasis achievement in humans



*When life threatening bleeding occurs, use QuikClot Combat Gauze® LE to stop bleeding quickly, efficiently, and effectively.*

## Unique Differentiator: Kaolin

*Kaolin is the product differentiator for all QuikClot products.*

QuikClot products...

- Use the body's natural clotting ability<sup>1</sup>
- Absorb blood and initiate the pathway to clot formation, leading to faster bleeding control<sup>1,3</sup>
- Include indications for external use to control traumatic bleeding
- Have consistent data to support product efficacy<sup>4,5,6</sup>

## Appropriate for traumatic bleeding injuries

- **Controls traumatic bleeding fast at the scene.**  
Promotes clotting within minutes, allowing you to control bleeding in an efficient and timely manner.<sup>3</sup>
- **Easy to use**  
Intuitive to use, QuikClot uses the same placement technique as standard gauze.
- **Updated packaging, easier to open**  
Compact, cube-sized packaging, and eight tear points expedite access to dressing for quicker application.

# QUICK. EFFICIENT. EFFECTIVE.

# Ordering Information

Product	Order No.	Description
QuikClot Combat Gauze® LE	350	Z-Fold (3 inch x 4 yards)
QuikClot® Belt Trauma Kit® (BTK)	261	1 Combat Gauze LE, 1 SWAT-T™, 1 Pair of Gloves
QuikClot® Belt Trauma Kit® (BTK) w/orange SOFTT-W™	534	1 Combat Gauze LE, 1 Mini Sheer, 1 Pair of Gloves, 1 Orange SOFTT-W™
QuikClot Bleeding Control Kit (BCK)	539	3 Combat Gauze LE, 3 Tourniquets, 1 Mini Sheer, 3 Mini Dressing
QuikClot Combat Gauze® LE Moulage Trainer	273	Z-Fold (3 inch x 4 yards)
QuikClot® Hemorrhage Control Training Kit (with QuikClot Combat Gauze® LE)	384	Visit QuikClot.com/Law-Enforcement for product details
QuikClot® Hemorrhage Control Training Kit with BioFeedback	529	Visit QuikClot.com/Law-Enforcement for product details

**Questions? Please contact us by Phone Mon-Fri 8am-5pm [ET]  
203.294.0000 or Email at [contactus@Z-medica.com](mailto:contactus@Z-medica.com)  
For additional information, visit us at [QuikClot.com/Law-Enforcement](http://QuikClot.com/Law-Enforcement)**

## References:

1. Lamb KM, Pitcher HT, Cavarocchi NC, Hirose H. Vascular Site Hemostasis in Percutaneous Extracorporeal Membrane Oxygenation Therapy. *Open Cardiovasc Thorac Surg J.* 2012;5:8-10
2. Trabattoni D, Montorsi P, Fabbicocchi F, Lualdi A, Gatto P, Bartorelli AL. A new kaolin-based haemostatic bandage compared with manual compression for bleeding control after percutaneous coronary procedures. *Eur Radiol.* 2011;21:1687-1691
3. Garcia-Blanco J, Gegel B, Burgert J, Johnson S, Johnson D. The Effects of Movement on Hemorrhage When QuikClot Combat Gauze® Is Used in a Hypothermic Hemodiluted Porcine Model. *J Spec Oper Med.* 2015 Spring;15(1):57-60
4. Zietlow JM, Zietlow SP, Morris DS, Berns KS, Jenkins DH. Prehospital Use of Hemostatic Bandages and Tourniquets: Translation From Military Experience to Implementation in Civilian Trauma Care. *J Spec Oper Med.* 2015 Summer;15(2):48-53.
5. Schauer SG, April MD, Naylor JF, Fisher AD, Cunningham CW, Ryan KL, Thomas KC, Brillhart DB, Fernandez JRD, Antonacci MA. QuikClot Combat Gauze® Use by Ground Forces in Afghanistan The Prehospital Trauma Registry Experience. *J Spec Oper Med.* Summer 2017;17(2):101-106.
6. Shina AT, Lipsky AM, Nadler R, Levi M, Benov A, Ran Y, Yitzhak A, Glassberg E. Prehospital use of hemostatic dressings by the Israel Defense Forces Medical Corps: A case series of 122 patients. *J Trauma Acute Care Surg.* 2015 Oct;79(4 Suppl 2):S204-9.

## Z-MEDICA®

4 Fairfield Boulevard, Wallingford, CT 06492  
Z-Medica.com

© 2019 Z-MEDICA, LLC. All rights reserved.  
PSST-1012